

# OS30p<sub>+</sub> Research Chlorophyll Fluorometer

## *The best of both worlds*

Hand held pulse modulated fluorometer for  $F_V/F_M$  /  $F_V/F_O$ , along with an advanced **OJIP** protocol with auto-calibration of actinic light source at  $3,500 \mu\text{mols m}^{-2}\text{s}^{-1}$



### Pulse Modulated Tests:


$F_V/F_M$  - A Measure of Maximum Quantum Efficiency

$F_V/F_M$  is by far the most used chlorophyll fluorescence measuring parameter in the world. Research has proven it to be a robust way to measure early plant stress that affects photosystem II. The parameter and its protocol also provide the advantage that samples are measured and compared at the same known dark adapted state.

$F_V/F_M$  has been shown to correlate with carbon assimilation for many types of plant stress in  $C_3$  plants,  $C_4$  plants & CAM plants. The OS30p<sub>+</sub> measuring trace is graphically displayed in color, and  $F_O$  is accurately measured using red modulated light.

$F_V/F_O$  - While it does not directly correlate with carbon assimilation, it is a very sensitive stress detector that is more sensitive than  $F_V/F_M$ . It also allows comparison of samples in the same known dark adapted state.

### The “JIP” test, with calibrated actinic light source - **Advanced OJIP**



Enhanced Strasser OJIP

**“JIP” Test** - Direct readout of important stress detecting parameters and the overlay of measuring traces, are now immediately possible in the field. This plant stress testing method provides a high time resolution image of the Kautsky induction curve against a logarithmic time scale. It is possible to display and overlay curve traces for comparison. They may also be easily recreated from the data file information to evaluate plant stress.

The OS30p<sub>+</sub> provides a direct read out of the following parameters: OJIP, t100 :s, t300 :s (or K), tFm (or time to  $F_M$ ), A (or area above the curve),  $M_O$  (or RC/ABS),  $PI_{ABS}$  (or performance Index),  $F_O$ ,  $F_M$ ,  $F_V/F_M$ , and  $F_V/F_O$ . In addition, OptiSciences goes the extra mile, by measuring  $F_O$  instead of estimating the parameter.

## $F_V/F_M$ - The most accepted & used chlorophyll fluorescence parameter in the world.

The reason why it is the most accepted chlorophyll fluorescence measuring protocol is because results correlate to carbon assimilation measurements and because all samples are compared at the same known dark adapted state.

$F_V/F_M$  is a normalized ratio that requires a modulated fluorometer for reliable measurement. The protocol was developed by Kitajima and Butler (1975) and it has withstood the test of time for measuring maximum quantum efficiency of PSII in plants. It has shown the robust capability to measure plant stresses that affect PSII.

After dark adaptation, usually with dark clips, modulated fluorometers allow the accurate measurement of minimum fluorescence or  $F_O$ . This is done by using a weak modulated light, that is too low to drive photosynthesis, but high enough to excite pre-photosynthetic antenna fluorescence. In this state, photosystem II is maximally oxidized. The xanthophyll cycle,  $\Delta$  pH of the thylakoid lumen, state transitions and chloroplast migration have all relaxed to their inactive states (Lichtenthaler 1999, 2004), (OptiSciences dark adaptation application note available upon request at [www.optisci.com](http://www.optisci.com)).


The modulated fluorometer then irradiates the plant sample with an intense saturation light that is both high and long enough, to fully reduce or close all available PSII reaction centers. The maximum fluorescence output during this saturating light radiation called  $F_M$  represents a fully reduced PSII.

It has been found that healthy plants have an  $F_V/F_M$  value in the range of 0.79 to 0.84 (Maxwell and Johnson 2000). Lower values indicate plant stress.

To make sure that the OS30p+ is the easiest to use and most reliable instrument in its class, it has the following design advantages:

1. Modulated light intensity setup bar.
2. An algorithm that finds the highest  $F_M$  as long as the saturation flash is wide enough. This eliminates saturation pulse NPQ as an issue.
3. A calibrated saturation flash up to  $6,000 \text{ umols m}^{-2}\text{s}^{-2}$  for reliable results under all temperature conditions.

## Graph of $F_V/F_M$ or $F_V/F_O$ saturation pulse



$$F_V/F_M = (F_M - F_O) / F_M$$

$$F_V/F_O = (F_M - F_O) / F_O$$

**No more getting up way before dawn to make field measurements!**



Inexpensive, long lasting, very light weight dark adaption clips allow measurement of statistically significant plant populations very quickly and at reasonable times of day. Each measurement takes between 3-5 seconds. *Made of light weight plastic & stainless steel components.*

*Each instrument comes with 10 dark clips. However, more may be purchase as singles or in packages of 50.*



## The “JIP” Test - OJIP

OJIP or the “JIP” test is another dark adapted test that has been used for detecting and measuring plant stress. It was discovered by Kautsky (1957). He used high time resolution of the rise in variable chlorophyll fluorescence, after dark adaption, using a saturating actinic light. He found a distinct curve shape with multiple steps that represent chemical functions in PSII oxidation-reduction reactions. Using this approach, plant stress that affects PSII can be measured.



The latest research shows that O-J is caused by photochemical quenching, J-I is caused by photoelectrochemical quenching, and I-P is associated with the electric trans-thylakoid potential generated by the proton pump fueled by Cyclic Electron Transport (CET) in PSI (Vredenberg 2011).

It has also been shown that some types of plant stress affect specific parts of the OJIP curve. In addition, special measuring parameters have been developed as sensitive stress detectors such as  $PI_{ABS}$  or performance index.

Viewing graphic OJIP results can now be quickly and easily done in the field. The OS30p<sub>+</sub> provides a color graphic display of the OJIP curve with a logarithmic time scale. It is common for researchers that use this technique to overlay measuring graph traces to study the effects of plant stress, and to use the special parameters that have been created to quantify plant stress. Up to 16 traces may be overlaid on the graphic color instrument screen, and up to 32 can be overlaid from a single measuring file in software.

Reto Strasser did much of pioneering work using OJIP for plant stress measurement.

Importantly, Wim Vredenberg (2011) found that OJIP results vary by actinic light intensity. Strasser originally used  $3,000 \mu\text{mols m}^{-2} \text{s}^{-1}$  in his early work, but he used  $3,500 \mu\text{mols m}^{-2} \text{s}^{-1}$  in his later work. For that reason, The OS30p<sub>+</sub> actinic light source goes through an auto-calibration routine every time measuring protocol selection occurs. This ensures reliable measurements at every temperature. The OS30p<sub>+</sub> may be used at either  $3,000 \mu\text{mols m}^{-2} \text{s}^{-1}$  or  $3,500 \mu\text{mols m}^{-2} \text{s}^{-1}$




The parameters O, J, I, P, t100us, t300us (or K),  $M_O$  (or RC/ABS),  $PI_{ABS}$  (or Performance Index), A (or Area above the curve), and tF<sub>M</sub> (or time to F<sub>M</sub>) are all displayed, along with  $F_v/F_m$ , and  $F_v/F_o$  on the screen.

Overlaid OJIP graphs and Spider graphs are now easily created using standard Microsoft Excel software. The data files are specifically organized with detailed time stamps, and consecutive row organization for very easy parameter spider graphing. Detailed OJIP overlay graphing is equally as easy, with the first column providing data capture time, and all OJIP curve traces lined up next to each other in columns to the right of the time. Both linear and logarithmic graphing become very easy to create.

USB output allows the OS30p<sub>+</sub>, to become a hard drive for your computer for simple transfer of data. No special external computer software is necessary.

The instrument measuring screen on the right allows instant review of the most important information when taking a measurement.


- Settings are shown in the upper left window
- A color graph of the measuring trace is shown using a logarithmic scale or a linear scale in the upper right window.
- Direct read out of the most used OJIP measuring parameters is shown in the lower window.



The last 100 measurements can be reviewed and compared as shown on screen to the right.


“JIP” test traces can now be view in the field without using a computer. Up to 32 traces can be overlaid on the viewing screen or in the data file.

Colors will start to repeat after 16 traces have been overlaid.




The default calibrated red light actinic intensity is 3,500  $\mu\text{mols}$ , other intensities are also available.

- 525  $\mu\text{mols}$
- 875  $\mu\text{mols}$
- 1,000  $\mu\text{mols}$
- 2,500  $\mu\text{mols}$
- 3,000  $\mu\text{mols}$
- 3,500  $\mu\text{mols}$
- 4,500  $\mu\text{mols}$
- 6,000  $\mu\text{mols}$



The graph on the right was made by using a measuring data file that was output by USB port to a spreadsheet. Sampling rate times are listed in column A of the spreadsheet and each successive trace reports to the following columns, one after the other.

The X axis is in logarithmic time and the Y axis is in fluorescence units. Of course, the Y axis may also be normalized over  $F_O$ .





**Rugged field instrument designed for hand held use.**

The OS30p+ has been one of most cited research fluorometers in the world in recent years.



The system comes with 10 light weight, quality dark adaptation clips that are very affordable. Pricing of dark clips allows the purchase of large quantities to fit most application requirements & budgets.



Standard accessories included are a hard shell carrying case, a battery charger, 10 dark adaptation clips, a manual on a USB drive and a USB cable.

## OS30p+ user references

### Cold stress:

Pokluda R., Sękara A., Jezdinsky A., Kalisz A. (2016) The physiological status and stress biomarker concentration of *Coriandrum sativum* L. plants subjected to chilling are modified by biostimulant application. *Biological Agriculture & Horticulture*: Volume 32, 2016 - Issue 4  
<https://www.tandfonline.com/doi/abs/10.1080/01448765.2016.1172344>

### Affects of LED growth lights:

Metallo R.M., Kopsell D.A., Sams C.E., Bumgarner N.R. (2018) Influence of blue/red vs. white LED light treatments on biomass, shoot morphology, and quality parameters of hydroponically grown kale. *Scientia Horticulturae*, Volume 235, 17 May 2018, Pages 189-197  
<https://www.sciencedirect.com/science/article/pii/S0304423818301456>

### Light stress under water deficit in papaya:

Vincent C., Schaffer B., Rowland D. (2018) Water-deficit priming of papaya reduces high-light stress through oxidation avoidance rather than anti-oxidant activity. *Environmental and Experimental Botany* Volume 156, December 2018, Pages 106-119  
<https://www.sciencedirect.com/science/article/abs/pii/S0098847218300893>

### Using urban organic waste in soil:

López-Cano I., Cayuela M., Sánchez-García M. (2018) —Suitability of different agricultural and urban organic wastes as feedstocks for the production of biochar Part 2: Agronomical evaluation as soil amendment *Sustainability*, 2018 - mdpi.com Volume 10, Issue 6 <https://www.mdpi.com/2071-1050/10/6/2077>

### The affects of heat and drought on photoprotective and antioxidant responses:

Wujeska-Klaus A., Bossinger G., Tausz M. (2015) Seedlings of two *Acacia* species from contrasting habitats show different photoprotective and antioxidative responses to drought and heatwaves. *Annals of Forest Science*. - Springer June 2015, Volume 72, Issue 4, pp 403–414

## OS30p+ user references

### The affects of microbial endophytes on conifers under stress:

MM Aghai, Z Khan, AM Stoda, AW Sher (2019) The Effect of Microbial Endophyte Consortia on *Pseudotsuga menziesii* and *Thuja plicata* Survival, Growth, and Physiology Across Edaphic Gradients. *Front. Microbiol.*, 18 June 2019 | <https://doi.org/10.3389/fmicb.2019.01353>

### The affects of saline and alkaline water on grass:

Cicchelli F.D.F, Wehr J.B., Dalzell S.A., Li C. (2016) Overhead-irrigation with saline and alkaline water: Deleterious effects on foliage of Rhodes grass and leucaena. *Agricultural Water Management* Volume 169, May 2016, Pages 173-182  
<https://www.sciencedirect.com/science/article/pii/S0378377416300786>

### The affect of cold stress on juvenile basil:

Kalisz A., Jezdinsky A., Pokluda R., Sękara A. (2016) Impacts of chilling on photosynthesis and chlorophyll pigment content in juvenile basil cultivars. *Horticulture, Environment, and Biotechnology* August 2016, Volume 57, Issue 4, pp 330–339  
K Thomason, MA Babar, JE Erickson, M Mulvaney (2018) Comparative physiological and metabolomics analysis of wheat (*Triticum aestivum* L.) following post-anthesis heat stress. *PLoS one*, 2018 - journals.plos.org Published: June 13, 2018  
<https://doi.org/10.1371/journal.pone.0197919>

### Using nano-iron III oxide particles to enhance germination, seedling growth, and salt tolerance in sorghum:

HF Maswada, M Djanaguiraman (2018) Seed treatment with nano-iron (III) oxide enhances germination, seedling growth and salinity tolerance of sorghum. *Journal of Agronomy and Crop Science* Volume 204, Issue 6  
<https://onlinelibrary.wiley.com/doi/abs/10.1111/jac.12280>

### Genetic response to cold stress in oil seed plants

D Horvath, JV Anderson, WS Chao, P Zheng (2019) Genes associated with chloroplasts and hormone-signaling, and transcription factors other than CBFs are associated with differential survival after low temperature treatments of *Camelina sativa* biotypes *PLoS one*, 2019 - journals.plos.org Published: May 31, 2019  
<https://doi.org/10.1371/journal.pone.0217692>



## OS30p+ user references

### Tomato - lettuce intercropping:

TPL da Cunha-Chiamolera (2017) Gas exchange, photosynthetic pigments, and growth in tomato: lettuce intercropping. Chilean journal of agricultural research versión On-line ISSN 0718-5839 Chil. j. agric. res. vol.77 no.4 Chillán 2017 <http://dx.doi.org/10.4067/S0718-58392017000400295>

### profiling of herbicide target-site genes in tumbleweed:

VK Varanasi, S Bayramov, PVV Prasad (2017) Expression Profiles of psbA, ALS, EPSPS, and Other Chloroplastic Genes in Response to PSII-, ALS-, and EPSPS-Inhibitor Treatments in *Kochia scoparia*. American Journal of ..., 2017 - researchgate.net

### Using magnetically treated water on Tabasco pepper, radish and corn:

DI Ospina-Salazar, JA Benavides-Bolaños (2018) Photosynthesis and biomass yield in Tabasco pepper, radish and maize subjected to magnetically treated water. *Corpoica Ciencia y ...*, 2018 - scielo.org.co

### Herbaceous plant ageing:

BM Quarles, DA Roach (2019) Ageing in an herbaceous plant: Increases in mortality and decreases in physiology and seed mass. *Journal of Ecology*, 2019 - Wiley Online Library

### Using sorgaab to control weeds in sesame:

JC Murimwa, JT Rugare, S Mabasa (2019) Allelopathic Effects of Aqueous Extracts of Sorghum (*Sorghum bicolor* L. Moench) on the Early Seedling Growth of Sesame (*Sesamum indicum* L.) Varieties and... *International Journal of Agronomy*, Volume 2019, Article ID 5494756, 12 pages <https://doi.org/10.1155/2019/5494756>

### The effects of salt stress on *Medicago truncatula*:

R Najar, S Aydi, S Sassi-Aydi, A Zarai (2019) Effect of salt stress on photosynthesis and chlorophyll fluorescence in *Medicago truncatula*. *Plant Biosystems-An*, 2019 - Taylor & Francis

## OS30p+ user references

### Heat stress in Eucalyptus trees:

A Wujeska-Klause, G Bossinger, M Tausz (2015) Responses to heatwaves of gas exchange, chlorophyll fluorescence and antioxidants ascorbic acid and glutathione in congeneric pairs of *Acacia* and *Eucalyptus*. *Trees* December 2015, Volume 29, Issue6, pp 1929–1941

### Used with lidar measurements:

AJ Maguire, JUH Eitel, LA Vierling, DM Johnson (2019) Terrestrial lidar scanning reveals fine-scale linkages between microstructure and photosynthetic functioning of small-stature spruce trees at the A Jezdinský, BP Konvičná (2019) forest-tundra ecotone. *Agricultural and Forest Meteorology –Volumes 269270*, 15 May 2019, Pages 157-168 Sledování intenzity fotosyntézy u skleníkové kultury zeleniny - theses.cz

### Drought stress in a C4 plant.

DB Medeiros, EC Silva, RJMC Nogueira (2013) Physiological limitations in two sugarcane varieties under water suppression and after recovering. *Theoretical and Experimental Plant Physiology On-line version ISSN 2197-0025 Theor. Exp. Plant Physiol.* vol.25 no.3 Campo dos Goytacazes 2013 Utilizing Ground Level Remote Sensing to Monitor Peatland Disturbance CN McCann - 2016 - macsphere.mcmaster.ca

### Drought stress on Soybean:

Vu Ngoc Thang, Tran Anh Tuan and Vu Dinh Hoa (2008) of water stress during soybean reproductive stages under nethouse conditions. *Journal of Science and Development 2008: Volume VI, Số 2: 116-121 HANOI UNIVERSITY OF AGRICULTURE*

### The effects of brushing golf green turf grass:

Gu, Chenchen (2016) The effect of brushing on creeping bentgrass putting green quality. [http://rave.ohiolink.edu/etdc/view?acc\\_num=osu1460571821](http://rave.ohiolink.edu/etdc/view?acc_num=osu1460571821)

### Silicon mitigation of boron deficiency in cotton:

Jonas Pereira de Souza Junior, Renato de Mello Prado, Marcilene Machado dos Santos Sarah Guilherme Felisberto (2019) Silicon mitigates boron deficiency and toxicity in cotton cultivated in nutrient solution. First published: 09 July 2019 <https://doi.org/10.1002/jpln.201800398>

## Technical Specifications:

### OS30p +

### Modulated Fluorometer

#### $F_V/F_M$ , $F_V/F_O$

Saturation intensity 600- 6000  $\mu\text{mols}$  set from 10% to 100%  
Saturation light source Array of red LEDs 660 nm.  
Modulated light source red 0.2 to 1.0  $\mu\text{mols}$ .

Detection method Pulse modulated

Detector and Filter Pin photodiode with 700-750 nm b/p filter

$F_V/F_M$  Test duration 0.1 seconds to 1.5 seconds.  
The default saturation pulse duration is set at 1.0 second;  
however, the software takes a rolling eight point 25 ms average  
to determine  $F_O$ , and  $F_M$ , making it ideal for both algae and land  
plants.

Modulated light manual adjustment from 30% to 80%

Parameters measured and graphed  $F_O$ ,  $F_M$ ,  $F_V/F_M$ ,  $F_V/F_O$

### JIP Test

Actinic light intensities - 6000  $\mu\text{mols}$ , 4500  $\mu\text{mol}$ , **3500  $\mu\text{mols}$** ,  
**3000  $\mu\text{mols}$** , 2500  $\mu\text{mols}$ , 1000  $\mu\text{mols}$ , 875  $\mu\text{mols}$ , 525  $\mu\text{mols}$ .  
An array of red LEDs at 650 nm are used for actinic illumination.

Detection method - Pin photodiode detector with 700-750 nm  
band pass filter using red pulse modulated light, and variable  
sampling rate from 10  $\mu\text{s}$  to seconds.

Test duration - JIP test 3 - 300 seconds

"JIP " Test parameters measured - Direct readout parameters:  
 $O$ ,  $t_{100\mu\text{s}}$ ,  $t_{300\mu\text{s}}$  (or K),  $t_{2\text{ms}}$  (or J),  $t_{30\text{ms}}$  (or I), P,  $t_{F_M}$ , A  
(area above the curve),  $M_O$  (or RC/ABS),  $PI_{ABS}$  (or performance  
index)  $F_O$ ,  $F_M$ ,  $F_V/F_M$ ,  $F_V/F_O$ ,

Up to 32 OJIP traces can be overlaid and displayed on the  
graphic display screen, Colors on the screen start to repeat after  
16 traces.  $F_o$  is measured, not estimated.

Measured parameters reported to the data file: ABS/RC,  
 $TR_O/RC$ ,  $DI_O/CS$ ,  $ET_O/RC$ ,  $TR_O/ABS$ ,  $ET_O/TR_O$ ,  $ET_O/CS$ ,  
 $RC/CS_O$ ,  $RC/CS_M$ , S, M, T are also measured and recorded in  
the data file, but not on the measuring screen.

Up to 32 traces can be overlaid on the graphic display screen,  
and up to 32 traces, per file name, can be stored for graphic  
overlay of traces in a single data file. Colors start to repeat on  
the viewing screen after 16 traces are overlaid. This limitation  
does not exist in the data file. If traces are not stored, thousands  
of measurement parameter sets can be stored in a single data  
file for spider graphing. The number of data files are only limited  
by machine memory limits.

## General Specifications

Auto Calibration of actinic light source when opening measuring  
protocols: Ensures that the actinic light intensity is the same even under  
hot or cold measuring conditions. It also ensures that OJIP  
measurements all occur at referenced actinic intensities.

Display - Color Graphic Display

Storage Capacity - Up to 160,000 measurements. Up to 32 traces can  
be stored with a one data file name. Hundreds of traces can be stored in  
multiple data files, with different names.

Digital output - USB port

Battery - NiMH battery pack with a battery life of 8 hours between  
charges.

Dimensions - 18cm, 7 cm, 6cm.

Weight - 1.25 lbs. , with carrying case and accessories - 4 lbs.

Carrying case - Included std. Case dimensions: 36cm, 28cm,  
15 cm.

### Technical references:

Vredenberg Wim (2011) Kinetic analyses and mathematical modeling  
of primary photochemical and photoelectrochemical processes in plant  
photosystems, BioSystems Contents lists available at ScienceDirect  
journal homepage: [www.elsevier.com/locate/biosystems](http://www.elsevier.com/locate/biosystems)

Strasser R.J, Tsimilli-Michael M., and Srivastava A. (2004) - Analysis  
of Chlorophyll a Fluorescence Transient. From Chapter 12,  
"Chlorophyll a Fluorescence a Signature of Photosynthesis", edited by  
George Papaqeorgiou and Govindjee, published by Springer 2004, PO  
Box 17, 3300 AA Dordrecht, The Netherlands, page 340

Kautsky H., Hirsch A. (1931) Neuw Versuche zur  
Kohlensaureassimilation. Naturwissenschaften 19, 964.

Baker N. R., Oxborough K., (2004) Chlorophyll fluorescence as a  
probe of photosynthetic productivity. From Chapter 3, "Chlorophyll a  
Fluorescence a Signature of Photosynthesis", edited by George  
Papaqeorgiou and Govindjee, published by Springer 2004, PO Box 17,  
3300 AA Dordrecht, The Netherlands, pages 66-79

Kitajima M, Butler WL (1975) Quenching of chlorophyll fluorescence  
and primary photochemistry in chloroplasts by dibromothymoquinone.  
Biochim Biophys Acta 376:105-115

Lichtenthaler H. K., Burkart S., (1999) Photosynthesis and high light  
stress. Bulg. J. Plant Physiol., 1999, 25(3-4), 3-16

Lichtenthaler H. K., Babani F. (2004) Light Adaption and Senescence  
of the Photosynthetic Apparatus. Changes in Pigment Composition,  
Chlorophyll Fluorescence Parameters and Photosynthetic Activity.  
From Chapter 28, "Chlorophyll a Fluorescence a Signature of  
Photosynthesis", edited by George Papaqeorgiou and Govindjee,  
published by Springer 2004, PO Box 17, 3300 AA Dordrecht, The  
Netherlands, page 716

Maxwell K., Johnson G. N, (2000) Chlorophyll fluorescence – a  
practical guide. Journal of Experimental Botany Vol. 51, No. 345, pp.  
659-668- April 2000